



Resuspension Reagent

Revision A, January 2024

Product No. KIT230010

Reagent for 96 tests

Store the kit at 15 to 25 °C

For food testing purposes.

FOR *IN VITRO* USE ONLY



Table of Contents

1. What this Product Does	3
1.1 Contents.....	3
1.2 Storage and Stability	3
2. How to Use this Product	3
3. Warranty and Disclaimer of Liability.....	3
4. Supplementary Information	4
4.1 Ordering Information.....	4
4.2 Trademarks.....	4
4.3 Contact and Support.....	4
4.4 Reference Number.....	4
5. Change Index	4



1. What this Product Does

The Resuspension Reagent is designed for the suspension of microbes in combination with the foodproof StarPrep Two Kit. The Reagent eliminates the need for hazardous organic extractions or chaotropic agents.

1.1 Contents

Product	Content	Storage
KIT230010	1 container with 35 mL Resuspension Reagent	15 to 25 °C

1.2 Storage and Stability

Store the Resuspension Reagent at room temperature (15 - 25 °C) through the expiration date printed on the label.

2. How to Use this Product

The Resuspension Reagent should be used in combination with the foodproof® StarPrep Two Kit (Product No. KIT230177) to provide the rapid preparation of bacterial DNA for direct use in PCR.

3. Warranty and Disclaimer of Liability

“Limited Warranty” and “Disclaimer of Liability”: Hygiena® Diagnostics warrants that this product is free from defects in materials and workmanship through the expiration date printed on the label and only if the following are complied with:

- (1) The product is used according to the guidelines and instructions set forth in the product literature;
- (2) Hygiena Diagnostics does not warrant its product against any and all defects when: the defect is as a result of material or workmanship not provided by Hygiena Diagnostics; defects caused by misuse or use contrary to the instructions supplied, or if the product is contaminated by improper storage or handling;
- (3) All warranties of merchantability and fitness for a particular purpose, written, oral, expressed or implied, shall extend only for a period of one year from the date of manufacture. There are no other warranties that extend beyond those described on the face of this warranty;
- (4) Hygiena Diagnostics does not undertake responsibility to any purchaser of its product for any undertaking, representation or warranty made by any dealers or distributors selling its products beyond those herein expressly expressed unless expressed in writing by an officer of Hygiena Diagnostics;
- (5) Hygiena Diagnostics does not assume responsibility for incidental or consequential damages, including, but not limited to responsibility for loss of use of this product, removal or replacement labor, loss of time, inconvenience, expense for telephone calls, shipping expenses, loss or damage to property or loss of revenue, personal injuries or wrongful death;
- (6) Hygiena Diagnostics reserves the right to replace or allow credit for any modules returned under this warranty.



4. Supplementary Information

4.1 Ordering Information

Hygiena Diagnostics offers a broad range of reagents and services. For a complete overview and for more information, please visit our website at www.hygiena.com.

4.2 Trademarks

foodproof® is a registered trademark of Hygiena Diagnostics GmbH. Other brand or product names are trademarks of their respective holders.

4.3 Contact and Support

If you have questions or experience problems with this or any other product of Hygiena Diagnostics GmbH, please contact our Technical Support staff (www.hygiena.com/support). Our scientists commit themselves to providing rapid and effective help. We also want you to contact us if you have suggestions for enhancing our product performance or using our products in new or specialized ways. Such customer information has repeatedly proven invaluable to us and the worldwide research community.

4.4 Reference Number

The reference number and original Hygiena Diagnostics GmbH article number: A 500 14

5. Change Index

Version 1, July 2014

First version of the package insert.

Revision A, January 2024

Rebranding and new layout.

A 500 14 20-> INS-KIT230010-RevA



Hygiena[®]

Camarillo, CA 93012

USA

diagnostics.support@hygiena.com

Manufactured by

Hygiena Diagnostics GmbH

Hermannswerder 17

14473 Potsdam

Germany

www.hygiena.com