




Stick Sponge™ | Sample Collection

Hygiena's Stick Sponges are biocide-free polyurethane sponges attached to a plastic handle and supplied in a sterile, tie-off bag. The easy-to-release stick makes it easier to reach into drains, pipes, and around equipment.

Benefits

- For surface sampling and hard to reach areas
- 1.5" x 3" biocide-free polyurethane sponge
- Sterilized by gamma irradiation
- Pre-moistened with a variety of diluents
- Sample collection without directly touching the sponge
- Easy to open tie-off bag
- Large write-on area for easy sample identification
- Attached sterile gloves optional

Stick Sponge Procedure



Product No.	Description	Quantity
SS100BPWP	Stick Sponge, 1.5" x 3" Sponge w/10 mL Buffered Peptone Water Broth, 5 bags of 20	100
SS100BPWPG	Stick Sponge, 1.5" x 3" Sponge w/10 mL Buffered Peptone Water Broth, 5 bags of 20 with gloves	100
SS100NBP	Stick Sponge, 1.5" x 3" Sponge w/10 mL Neutralizing Buffer, 5 bags of 20	100
SS100NBPG	Stick Sponge, 1.5" x 3" Sponge w/10 mL Neutralizing Buffer, 5 bags of 20 with gloves	100
SS100LETP	Stick Sponge, 1.5" x 3" Sponge w/10 mL Lethen Broth, 5 bags of 20	100
SS100LETPG	Stick Sponge, 1.5" x 3" Sponge w/10 mL Lethen Broth, 5 bags of 20 with gloves	100
SS100DEP	Stick Sponge, 1.5" x 3" Polyurethane Sponge w/10mL D/E Neutralizing Broth, 24 oz., 5 bags of 20	100
SS100DEPG	Stick Sponge, 1.5" x 3" Polyurethane Sponge w/10mL D/E Neutralizing Broth, 24 oz., 5 bags of 20 with gloves	100




Sponge'n Bag™ | Sample Collection

Hygiena's Sponge'n Bags are biocide-free polyurethane sponges supplied in a sterile, tie-off bag. The sponge provides a simple and convenient device for collecting environmental samples from large surface areas. Sponges may be ordered pre-moistened with a variety of different diluents.

Benefits

- 1.5" x 3" biocide-free polyurethane sponge
- Sterilized by gamma irradiation
- Pre-moistened with a variety of diluents
- Easy to open tie-off bag
- Large write-on area for easy sample identification
- Attached sterile gloves optional

Sponge'n Bag Procedure



Product No.	Description	Quantity
HSB100BPWP	Sponge'n Bag, 1.5" x 3" Sponge w/10 mL Buffered Peptone Water Broth, 5 bags of 20	100
HSB100BPWPG	Sponge'n Bag, 1.5" x 3" Sponge w/10 mL Buffered Peptone Water Broth, 5 bags of 20 with gloves	100
HSB100LETP	Sponge'n Bag, 1.5" x 3" Sponge w/10 mL Lethen Broth, 5 bags of 20	100
HSB100LETPG	Sponge'n Bag, 1.5" x 3" Sponge w/10 mL Lethen Broth, 5 bags of 20 with gloves	100
HSB100NBP	Sponge'n Bag, 1.5" x 3" Sponge w/10 mL Neutralizing Buffer, 5 bags of 20	100
HSB100NBPG	Sponge'n Bag, 1.5" x 3" Sponge w/10 mL Neutralizing Buffer, 5 bags of 20 with gloves	100
HSB100DEP	Sponge'n Bag, 1.5" x 3" Sponge w/10 mL D/E Neutralizing Broth, 5 bags of 20	100
HSB100DEPG	Sponge'n Bag, 1.5" x 3" Sponge w/10 mL D/E Neutralizing Broth, 5 bags of 20 with gloves	100